



# Chair *Caddie*

RETRO-FIT



***Turns...Rolls...and Brakes for Safety!***



Enables a standing person get seated... without the fear of falling.



Enables a caregiver to move a seated person up-to the table...with ease, grace, and safety!

- Extends life of chairs by 4-6 years
- Eliminates damage to floors
- Attaches to existing chairs



 [www.comfortek.com/cc](http://www.comfortek.com/cc)

- Reduces caregiver fatigue
- Reduces workplace injuries
- Provides dignity to persons with limited mobility

# Chair *Caddie*

RETRO-FIT

1. **Attaches to existing chairs** (wood or metal)
2. Low profile design – chair is raised  $\frac{3}{4}$ " above the floor
3. Fits inside the footprint of an existing chair
4. Medical grade casters:
  - Chair Caddie A (3" front / 2" rear)
  - Chair Caddie B (4" front / 3" rear)
5. Anti-tip glides
6. Front wheel tension control
7. Foot activated locking rear caster
8. Rated for 450 lbs.
9. Attaches to arm or armless chairs
10. Adjustable - Custom sizing available
11. Eliminates **pushing, pulling, shoving**, and **twisting**



PATENT  
PENDING



The Chair CADDIE's mobility features  
(**turn...roll...and brake** for safety) help to;

- **reduce** a standing person's risk of falling
- **reduce** caregiver's risk of incurring a work / personal injury!








APRIL 2023

The Chair CADDIE is an aftermarket accessory which attaches to the base of many of today's dining and banquet chairs as well as lobby/guest/side chairs!

By design, the Chair CADDIE is an "undercarriage unit" that provides mobility, structural integrity and stability to chairs, without compromising safety for both the caregiver or the seated person.


When attached, the chair along with the seated person **swivels...turns...rolls...**and **brakes** for safety, eliminating the need for caregiver to push-pull-shove-twist each time a seated person is moved up-to or away-from the table.

Operationally, the Chair CADDIE reduces the physical demands of the caregiver providing mealtime assistance in three (3) specific caregiving situations:



**1. Assist a person into the chair**


Roll the Chair CADDIE directly behind the person until the front of the chair seat touches the back of their legs. After the brakes are engaged, assist the person into the chair.



**2. Assist the seated person up-to the table**

With the person comfortably seated, the caregiver is able to gracefully roll the seated person up-to the table. For safety, the brakes can be applied once the person is seated at the table

**NOTE:** The walker/rollator should be moved off to the side of the dining room freeing the aisle for use by others.



**3. Assist the seated person away-from the table**

To exit the table, the seated person is effortlessly rolled away from the table and turned in the direction of the walker/rollator.

Once the walker/rollator has been retrieved and placed in front of the seated person, the seated person can be assisted to a standing position by the caregiver. Holding onto the walker/rollator for support the person walks directly away from the table.

## DESIGN ELEMENTS OF THE CHAIR CADDIE



1. The Chair CADDIE is a fully integrated “undercarriage unit” which strengthens and supports the frame of the chair to which it is attached.
2. Side rails hold the legs of the chair in place preventing them from splaying outwards. The connecting arms are adjustable enabling the Chair CADDIE to attach to most existing armed and armless chairs. The Chair CADDIE is rated to carry 450 lbs.
3. Casters (oversized) attached directly to the “undercarriage unit” (not to the legs of the chair) enable the chair with the person seated in the chair to move easily on tile, laminate or carpeted surfaces, while raising the seat height of the chair only 3/4”.
4. For stability, the front casters are positioned in-line with the front legs of the chair. The wide wheel-base stance provides unparalleled stability preventing the chair from tipping sideways. Anti-tip glides installed under the front legs prevent the chair from tipping forwards.
5. For mobility, the “undercarriage unit” is fitted with rotating rear casters.
6. For safety, both rear casters are fitted with foot activated brakes. It is recommended the brakes are engaged whenever the chair is not attended by a caregiver.
7. Additional, the entire “undercarriage unit” fits inside the “foot print” of the chair to which it is attached and does not pose a tripping hazard to the caregiver or other guests.



Online Reference: [www.comfortek.com/cc](http://www.comfortek.com/cc)

## OPERATING GUIDELINES



info@comfortek.com  
1.888.678.2060