

## T2-841-22-CC2

Action: Designed to turn/roll...LOCKABLE by User

Benefits: Designed for persons who live independently, but

experience difficulty "scooting" themselves upto the table. Persons can confidently be seated knowing the brakes prevent the chair from rolling away. Once seated, the seated person releases the brakes, enabling them to propel themselves (sit-walking) up-to the table using both hands/feet

for guidance.

Features: An armed chair (22" wide) to which the CADDIE is

attached. The Brake-Release-Lever is positioned beside the seat of the chair...in easy reach of the seated person. The CADDIE includes 4 casters, 2

of which rotate 360°.