



## *Chair Caddie A Retrofit*

Action: **Designed to...roll / LOCK (DIY)**

Benefits: The CADDIE enables a caregiver to move a seated person up-to and away-from the table without the need to push / pull / shove / twist.

Features: DIY...A metal platform complete with 4 inline casters which attaches to the base of an existing chair. One of the casters includes a foot / toe activated locking lever.

# Chair *Caddie*<sup>TM</sup>

RETRO-FIT



It Rolls, It Turns



Patent Application No. 15/996,068



Attaches to ***YOUR*** dining chairs



The CHAIR CADDIE enables a seated person to be moved up-to and away-from the table with **EASE** and **SAFETY!**

- ♥ Reduces caregiver fatigue / burnout
- ♥ Provides dignity to persons with limited mobility



Watch video at [www.comfortek.com/CC](http://www.comfortek.com/CC)

# Chair *Caddie*™

RETRO-FIT



1. **Attaches to existing chairs** (wood or metal)
2. Low profile design – chair is raised  $\frac{3}{4}$ " above the floor
3. Fits inside the footprint of an existing chair
4. Medical grade casters (3" / 4" Diameter)
5. Anti-tip glides
6. Front wheel tension control
7. Foot activated locking rear caster
8. Rated for 450 lbs.

**MODEL A** (Rolls...and LOCKS for safety)

(X1) 16" - 23"	(Y) 14 $\frac{1}{2}$ " - 23"	(Z) 15 $\frac{1}{2}$ " - 20 $\frac{1}{4}$ "
(X2) 13 $\frac{1}{8}$ " - 19 $\frac{3}{8}$ "		

**MODEL B** (Rolls...turns...and LOCKS for safety)

(X1) 17 $\frac{5}{8}$ " - 24"	(Y) 17" - 24"	(Z) 18" - 24 $\frac{1}{16}$ "
(X2) 11 $\frac{7}{16}$ " - 15 $\frac{5}{16}$ "		

Watch video at [www.comfortek.com/CCA](http://www.comfortek.com/CCA)

Watch video at [www.comfortek.com/CCB](http://www.comfortek.com/CCB)

Changing the *costs culture community* of care



www.comfortek.com  
1-888-678-2060