



## *Chair Caddie A Retrofit*

Action: **Designed to...roll / LOCK (DIY)**

Benefits: The CADDIE enables a caregiver to move a seated person up-to and away-from the table without the need to push / pull / shove / twist.

Features: DIY...A metal platform complete with 4 inline casters which attaches to the base of an existing chair. One of the casters includes a foot / toe activated locking lever.

# Chair *Caddie*<sup>TM</sup>

RETRO-FIT



It Rolls, It Turns



Attaches to ***YOUR*** dining chairs



The CHAIR CADDIE enables a seated person to be moved up-to and away-from the table with **EASE** and **SAFETY**!



*Reduces caregiver fatigue / burnout*



*Provides dignity to persons with limited mobility*



Watch video at [www.comfortek.com/CC](http://www.comfortek.com/CC)


# Chair *Caddie*™

RETRO-FIT



1. **Attaches to existing chairs** (wood or metal)
2. Low profile design – chair is raised  $\frac{3}{4}$ " above the floor
3. Fits inside the footprint of an existing chair
4. Medical grade casters (3" / 4" Diameter)
5. Anti-tip glides
6. Front wheel tension control
7. Foot activated locking rear caster
8. Rated for 450 lbs.

## MODEL A (Rolls...and LOCKS for safety)




(X1) 16" - 23"

(Y) 14  $\frac{1}{2}$ " - 23"

(Z) 15  $\frac{1}{2}$ " - 20  $\frac{1}{4}$ "

(X2) 13  $\frac{1}{8}$ " - 19  $\frac{3}{8}$ "

## MODEL B (Rolls...turns...and LOCKS for safety)



(X1) 17  $\frac{5}{8}$ " - 24"

(Y) 17" - 24"

(Z) 18" - 24  $\frac{1}{16}$ "

(X2) 11  $\frac{7}{16}$ " - 15  $\frac{5}{16}$ "



Watch video at [www.comfortek.com/CCA](http://www.comfortek.com/CCA)



Watch video at [www.comfortek.com/CCB](http://www.comfortek.com/CCB)

Changing the *costs culture community* of care



[www.comfortek.com](http://www.comfortek.com)  
1-888-678-2060