



Chair *Caddie*

RETRO-FIT



Rolls...and Brakes for Safety!



Enables a standing person get seated... without the fear of falling.



Enables a caregiver to move a seated person up-to the table...with ease, grace, and safety!

- Extends life of chairs by 4-6 years
- Eliminates damage to floors
- Attaches to existing chairs



- Reduces caregiver fatigue
- Reduces workplace injuries
- Provides dignity to persons with limited mobility

 www.comfortek.com/cc

Chair *Caddie*

RETRO-FIT

1. **Attaches to existing chairs** (wood or metal)
2. Low profile design – chair is raised $\frac{3}{4}$ " above the floor
3. Fits inside the footprint of an existing chair
4. Medical grade casters:
 - Chair Caddie A (3" front / 2" rear)
 - Chair Caddie B (4" front / 3" rear)
5. Anti-tip glides
6. Front wheel tension control
7. Foot activated locking rear caster
8. Rated for 450 lbs.
9. Attaches to arm or armless chairs
10. Adjustable - Custom sizing available
11. Eliminates **pushing, pulling, shoving, and twisting**



PATENT
PENDING



The mobility features of the Chair CADDIE (**roll...and brake** for safety) help to:

- **Prevent** resident falls
- **Reduce** care staff injuries
- **Increase** resident dignity and self-worth

