



Chair *Caddie*

RETRO-FIT



Turns...Rolls...and Brakes for Safety!



Enables a standing person get seated... without the fear of falling.



Enables a caregiver to move a seated person up-to the table...with ease, grace, and safety!

- Extends life of chairs by 4-6 years
- Eliminates damage to floors
- Attaches to existing chairs



 www.comfortek.com/cc

- Reduces caregiver fatigue
- Reduces workplace injuries
- Provides dignity to persons with limited mobility

Chair *Caddie*

RETRO-FIT

1. **Attaches to existing chairs** (wood or metal)
2. Low profile design – chair is raised $\frac{3}{4}$ " above the floor
3. Fits inside the footprint of an existing chair
4. Medical grade casters:
 - Chair Caddie A (3" front / 2" rear)
 - Chair Caddie B (4" front / 3" rear)
5. Anti-tip glides
6. Front wheel tension control
7. Foot activated locking rear caster
8. Rated for 450 lbs.
9. Attaches to arm or armless chairs
10. Adjustable - Custom sizing available
11. Eliminates **pushing, pulling, shoving, and twisting**



PATENT
PENDING



The Chair CADDIE's mobility features
(**turn...roll...and brake** for safety) help to;

- **reduce** a standing person's risk of falling
- **reduce** caregiver's risk of incurring a work / personal injury!

