



Chair **Caddie**TM

RETRO-FIT



It Turns...It Rolls...and Brakes for Safety!

- Extends life of chairs by 2-4 years
- Eliminates damage to floors
- Attaches to existing chairs



The Chair CADDIE enables a spouse / caregiver to easily and safely move a seated person up-to and away-from the table...without risk of injury.

- Reduces caregiver fatigue
- Reduces workplace injuries
- Provides dignity to persons with limited mobility

NO Pushing...NO Shoving...NO Twisting!



www.comfortek.com/cc

**PATENT
PENDING**

Chair Caddie™


RETRO-FIT

1. **Attaches to existing chairs** (wood or metal)
2. Low profile design – chair is raised $\frac{3}{4}$ " above the floor
3. Fits inside the footprint of an existing chair
4. Medical grade casters:
Chair Caddie A (3" front / 2" rear)
Chair Caddie B (4" front / 3" rear)
5. Anti-tip glides
6. Front wheel tension control
7. Foot activated locking rear caster
8. Rated for 450 lbs.
9. Attaches to arm or armless chairs
10. Adjustable - Custom sizes available



Chair Caddie A

- HomeCARE
- Residential



X1 15 $\frac{3}{8}$ " - 23 $\frac{5}{8}$ "


Y 14 $\frac{3}{4}$ " - 23"

Z 16 $\frac{7}{16}$ " - 22 $\frac{9}{16}$ "

X2 10 $\frac{11}{16}$ " - 18 $\frac{15}{16}$ "

Chair Caddie B

- Senior Living / LTC
- Adult Day Care
- Commercial



X1 17 $\frac{5}{8}$ " - 24"

Y 17" - 24"

Z 18" - 24 $\frac{1}{16}$ "

X2 13" - 19 $\frac{3}{8}$ "



Watch video at www.comfortek.com/cc

Changing the *costs* culture *of care* community



www.comfortek.com | 1-888-678-2060